

## *Eriochrome Black-T Indicator Color Changes*



Since the reactants and products in the titration between hard water and EDTA (equation 1) are all colorless, an indicator is needed to allow us to visually see the *end point* (or *equivalence point*) of the titration, i.e. *the point at which moles of  $\text{Ca}^{2+}$  = moles of EDTA*. Eriochrome Black-T is the indicator we will use.



### *How Eriochrome Black-T Works*

When Eriochrome Black-T complexes with  $\text{Mg}^{2+}$  ions, it produces a **PINK-RED** solution. The indicator is **BLUE** when it is not complexed with  $\text{Mg}^{2+}$  (equation 2) and the solution is basic.

As EDTA is added to hard water during the titration it reacts with the divalent ions (e.g.  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$ ) present in the hard water. Once all the divalent ions in solution have reacted with EDTA, the EDTA reacts with the  $\text{Mg}^{2+}$  ions that are complexed with the Eriochrome Black-T indicator, thus causing the solution to turn **BLUE**.




Click here for an [Animation of EDTA endpoint using Eriochrome Black-T Indicator](#)

The three pictures show the color changes as the end point of the titration is reached. There is only a 1 drop difference of 0.010 M EDTA between the first and second pictures and between the second and third pictures. Two or three seconds were allowed for colors in the second and third pictures to develop after adding the additional drop. In each case the solution was thoroughly mixed. This color change from wine red to violet to blue is due to the color changes of the indicator Eriochrome Black-T.

### Sample Problem

Suppose you are using **EDTA with a molarity of 0.008000** for the titration. You titrate **50.00 mL of a water sample** using **10.68 mL of 0.00800 M EDTA**. Calculate the concentration of  $\text{Ca}^{2+}$  ion in the water sample in *moles per liter* and in parts per million, *ppm  $\text{CaCO}_3$* .

**Answers: 0.001700 M  $\text{Ca}^{2+}$  or 170.0 ppm  $\text{CaCO}_3$**



$$\text{ppm} = \text{mass fraction} \times 10^6 \quad \text{or} \quad \text{ppm} = \left[ \frac{\text{mass}_{\text{solute}}}{\text{mass}_{\text{solution}}} \right] \times 10^6$$